

Introducing ProCair Sealed Base
by Novis – alternating solutions for
the prevention and treatment of pressure
injuries and other preventable complications
of immobility.

| | Standard size | King Single Size | |
|--|--|--|--|
| CODE | APMPC-R02S | APMPC-R02SK | |
| CAPACITY (SAFE WORKING LOAD) | 40 - 220 kg | | |
| SYSTEM ALTERNATION | 1 in 2 | | |
| WARRANTY | 3 year | | |
| CARRY BAG | Canvas Carry Bag | | |
| PRESSURE INJURY RISK | High to very high risk | | |
| COMPLIANCE | IEC 60601-1, IEC60601-1-2, IEC60601-1-11 | | |
| ARTG | 289458 | | |
| CARTON DIMENSIONS (L x W x H) in mm | 1120 x 430 x 400 mm | 1280 x 450 x 430 mm | |
| DIMENSIONS (L x W x H) in mm | 2000 x 880 x 200 mm | 2000 x 1050 x 200 mm | |
| MATTRESS WEIGHT | 13 kg | 14.5 kg | |
| WEIGHT SENSOR PAD | Yes | | |
| CPR TYPE | QuickTwist CPR at head end for rapid deflation | | |
| CELLS | NUMBER | 20 | |
| | STATIC HEAD CELLS / ZONE | 3 | |
| | NARROW FOOT CELLS / ZONE | 7 individually disconnectable | |
| | CELL-IN-CELL | Yes | |
| | CELL MATERIAL | Pure TPU | |
| | MATTRESS AIR SUPPLY LINE | Welded TPU air manifold and CPC quick release connectors | |
| TOP COVER | MATERIAL / TYPE | Dark blue, polyurethane-laminated nylon – water resistant, vapour permeable multi-stretch | |
| | FIRE RETARDANCY RATE | BS7175 Crib 5 | |
| | MAXIMUM WASH TEMPERATURE | 95° C | |
| | ATTACHEMENT | Dual zippers | |
| | SEAM CONSTRUCTION | RF welded on all external seams | |
| BASE | MATERIAL / TYPE | TPU laminated nylon | |
| | HINGE | No | |
| | FIXING STRAPS | Webbing at 4 corners, 2 at middle | |

Support
Surface

| | | |
|-------------------------------------|------------------|--|
| DIMENSIONS (L x W x H) in mm | | 350 x 135 x 224 mm |
| WEIGHT | | 3.4 kg |
| POWER | | 100-240V 50/60Hz |
| COMPRESSORS | | Single |
| OPERATION | | Full digital control |
| CYCLE TIME | | 12 minutes |
| NOISE LEVEL, WITHIN 2 METERS | | 32 db(A) ± 1 |
| MIN / MAX PRESSURE | | 20-60mm HG +/- 6mm HG |
| ELECTRICAL PROTECTION TYPE | | Class II Type BF |
| OPERATING ENVIRONMENT | | Operation 15°C to 35°C Storage 5°C to 60°C |
| PRESSURE SETTING | AUTOMATIC | Yes – constantly monitors weight and adjusts |
| | MANUAL | Low range, standard, high range |
| CAIRALERT PATIENT EXIT ALARM | | Yes |
| FOWLER FUNCTION | | Yes, automatic |
| BACKUP BATTERY | | Up to 5 hours continuous alternating therapy in the event of power failure. |
| ALARMS | | Low Pressure High Pressure Alternation Failure Startup Failure Patient Bottoming Out Sensor Disconnection Patient Exit Low Battery Service Due |

Control
Unit



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.